


DP-30Vet

Entry level system for veterinarians

BCF
Ultrasound



mindray
healthcare within reach

Ideal for Equine Reproduction

The entry level Mindray DP-30Vet is a lightweight, portable ultrasound system ideal for your equine or mixed practice. The Mindray DP-30Vet has been developed from the Mindray DP-50Vet and offers similar functionality including the same iClear speckle reduction technology.

Ergonomic design with fashionable outlook

- 12.1" LED monitor with tilting functionality
- Compact design for extreme mobility
- Fast boot-up within 17s for immediate scan
- User-defined keys, programmable personalized operations
- One-key function for image and cine storage to local disk or USB drive Lithium battery supporting up to 1.5 hours continuous scanning

Improved LED Screen

The new 12.1" LED Screen can be angled up or down, this is a really useful function when scanning at different heights.

iStation Integrated Patient Data Storage

iStation™ integrated patient data storage, archiving, review and retrieval efficiently.

Optional 320G hard disk for massive image storage

Specific design for versatile veterinary applications

- Dustproof keyboard cover keeps the control panel free from contamination
- Specifically designed veterinary transducers with anti-decay surface and longer cable

Transducers



65C15EA 5-8.5MHz
Convex Transducer



35C50EA 2-6MHz
Convex Transducer



35C20EA 2-5MHz
Convex Transducer



75L38EA 5-10MHz
Linear Transducer



75L50EAV 5-10MHz
Rectal Transducer



50L60EAV 4-6MHz
Rectal Transducer

Accessories



Padded Carry Bag



Battery Pack



Mobile Trolley



Power Lead

Specifications

System Weight	5.5kg without battery	Transducers	Range of transducers available
Dimensions	16.1cm x 29.0 cm x 35.4 cm	Colour Doppler	No
Display	12.1-inch LED High Resolution	Portable	Yes
Architecture	All digital broadband	Warranty	18 month warranty on system
Power	AC Power or lithium ion battery	Manufacturer	Mindray

Australia

t: 1300 798 747 | +61 3 9894 8980
sales@bcfultrasound.com

New Zealand

t: +64 9 257 1212
sales@bcfultrasound.com