


DP-50Vet

Ideal system for the general practitioner

BCFTM
Ultrasound

mindray
healthcare within reach



One system for all applications

Mindray's DP-50Vet black & white ultrasound system allows you to do more in less time without any compromise in quality of diagnosis. Designed as a stepping stone into the world of ultrasound, DP-50Vet will be your legend for general practices, while maintaining your budget needs.

Advanced imaging technologies

- iTouch™ Auto Optimization
- One touch image optimization, quickly optimizes the uniformity and dynamic range, very convenient and efficient

Exceptional ergonomics by design

- Sleek, streamlined, compact shape
- Lightweight system for enhanced mobility
- High resolution, wide-angle 15" LCD with tilt functionality for better viewing

iClear™ Speckle Suppression

Scale the image size to intelligently detect structures of different sizes, improve accuracy

Evaluating patients in today's demanding environment takes innovative technology to speed up diagnosis with greater confidence. Mindray's DP-50 versatility expands your clinical utility while improving patient care. The Mindray DP-50 is very easy to use and the quality of the ultrasound images produced is superior to several other more expensive machines which are available.

Transducers



65C15EA 5-8.5MHz
Convex Transducer



35C50EA 2-6MHz
Convex Transducer



35C20EA 2-5MHz
Convex Transducer



75L38EA 5-10MHz
Linear Transducer



75L50EAV 5-10MHz
Rectal Transducer



50L60EAV 4-6MHz
Rectal Transducer

Accessories



Padded Carry Bag



Battery Pack



Mobile Trolley



Standoffs

Specifications

System Weight	7.7kg without battery	Transducers	Range of transducers available
Dimensions	19 x 41.5 x 37.8cm	Colour Doppler	No
Display	15" diagonal LCD, 1024 x 768	Portable	Yes
Architecture	All digital broadband	Warranty	18 month warranty on system
Power	AC Power or lithium ion battery	Manufacturer	Mindray

Australia

t: 1300 798 747 | +61 3 9894 8980
sales@bcfultrasound.com

New Zealand

t: +64 9 257 1212
sales@bcfultrasound.com