


# M-Turbo

Focused on the Vet



Superior Image Quality



Robust Design

Portable and easy to use

## See More. Achieve More.

So tough, the VET M-Turbo system can withstand a 1 metre drop thanks to the unique magnesium frame. And for complete confidence every system comes with SonoSite's industry-leading 5-year warranty. Supremely reliable, easy to use, and durable enough to make it a favourite with vets everywhere – SonoSite sets the standard for field, barn, and exam room use.



## Applications



### See More. Do More.

Smaller than most laptops and weighing in at 3.04 kg, the VET M-Turbo hand-carried ultrasound system combines revolutionary image quality with portability in a robust casing.

## 5 Year Warranty - Industry First

For complete confidence every system comes with SonoSite's industry leading 5-year warranty. Supremely reliable, easy to use, and durable enough to make it a favourite with vets everywhere – SonoSite sets the standard for field, barn, and exam room use.

## Reliable

The M-Turbo has been designed to American military specifications and extensively tested to ensure that it is reliable and durable. So tough, the Vet M-Turbo system can withstand a 1 metre drop test thanks to the unique magnesium frame.

## Transducers - One for every application



L52x 10-5MHz 52mm  
Linear Rectal



HFL50x 15-6MHz 50mm  
Linear Array



L25x 13-6MHz 25mm  
Linear Array



HFL38x 13-6MHz 38mm  
Linear Array



C11x 8-5MHz 11mm  
Micro Convex Array



C60x 5-2MHz 60mm  
Convex Array



P10x 8-4MHz 10mm  
Phased Array



P21x 5-1MHz 21mm  
Phased Array



## Flexible

The M-Turbo is easily expandable, with an extensive range of transducers and functions including Colour, Doppler and Tissue Harmonic Imaging to suit any application.

## Excellent Image Quality

Using advanced software and hardware, the SonoSite M-Turbo provides revolutionary image quality.

## Portable


Similar in size and weight to a laptop computer, the SonoSite M-Turbo combines an elegant design with portability and a robust casing.

## Specifications

<b>System Weight</b>	3.04 kg with battery	<b>Architecture</b>	All digital broadband
<b>Dimensions</b>	30.2 cm x 27.4cm x 7.9 cm	<b>Dynamic Range</b>	Up to 165 dB
<b>Display</b>	10.4" / 26.4cm diagonal LCD	<b>Grey Scale</b>	256 Shades
<b>Image Storage</b>	8GB Internal Storage	<b>Power</b>	2hr Li Battery or AC Power
<b>Operation Modes</b>	B, B+B, B+M, M, CPD, PW, CW	<b>Grey Scale</b>	256 Shades
<b>Manufacturer</b>	Sonosite Fujifilm	<b>Warranty</b>	5 Year Warranty

## Optional Extras

BCF Ultrasound offer a range of accessories and software options for the M-Turbo Ultrasound System from gel standoffs to Continuous Wave Doppler



### Software

- Colour Power Doppler / Velocity Colour
- Continuous Wave Doppler
- Pulsed Wave Doppler
- SonoMBe
- Advanced OB/GYN Calculations
- Credentialing Worksheets
- DICOM
- ACO Technology
- SiteLink Image Manager 4.1.0
- Wireless Data Management
- ECG Option

### Accessories

Universal Stand	TTC
Mini-Dock	Carry Case
Power Supply	Extra Battery
L38/50 Standoff	L52 Standoff
Keypad Skins	Aquasonic Gel

## Training Courses

At BCF we like to know that once a product has left our warehouse, it is going to be used to its full potential. So it is our commitment to our customers to provide the training required to get the most out of your product. To achieve this, we hold and support CPD training courses throughout the year, as well as providing networks to assist continued education for all customers, hosting local and international talent to pass on their knowledge for the benefit of the industry. You can also check out our great range of online training videos.

## Service and Support

Whether you are after some useful information to help get the best out of your equipment or you are having some technical difficulties, we are here to assist you.

### Australia

t: 1300 798 747 | +61 3 9894 8980  
sales@bcfultrasound.com

### New Zealand

t: +64 9 257 1212  
sales@bcfultrasound.com